



## Mark 5 Plus Oxygen Concentrator Domestic Models –115V; 60Hz

Configurations 1373-31100 and 1373-71100 (Part 1)  
Configurations 1373-31103 and 1373-71103 (Part 2)

|  |                     |
|--|---------------------|
| Introduction.....                          | 1-1                 |
| Specifications.....                        | 1-1                 |
| Safety Information.....                    | 1-2                 |
| Principles of Operation.....               | 1-3                 |
| Pneumatic Diagram.....                     | 1-4                 |
| Component Function.....                    | 1-4                 |
| Initial Setup Instruction.....             | 1-5                 |
| Precautions and Safe Practices.....        | 1-6                 |
| Preventive Maintenance.....                | 1-7                 |
| Trouble Shooting.....                      | 1-8                 |
| Test Port Diagram.....                     | 1-8                 |
| System Diagnosis Chart.....                | 1-9                 |
| Trouble Shooting Guide.....                | 1-10                |
| Component Replacement.....                 | 1-11                |
| Reference Material                         |                     |
| Electrical Schematic.....                  | 1-13                |
| Illustrations and Part Numbers.....        | 1-14                |
| Oxygen Monitoring System Information       | ② See Section Two   |
| Price List/Sample Forms/Reference Material | ③ See Section Three |

Nidek Medical Products, Inc.  
3949 Valley East Industrial Drive  
Birmingham, AL 35217 U.S.A.  
(800) 822-9255  
Phone (205) 856-7200  
Fax (205) 856-0533





## Introduction

Thank you for purchasing the *Mark 5 Plus* Oxygen Concentrator.

The *Mark 5 Plus* Concentrator is the latest in the high quality Mark series of concentrators manufactured by Nidek Medical Products, Inc. It is a reliable concentrator, requiring very little maintenance and is designed to operate efficiently for a long period of time. We have accomplished this by combining real world experience with the latest of control and process technology.

We believe that the *Mark 5 Plus* is the best machine available in today's market and will operate to your complete satisfaction. We invite your comments on our devices and our service. The Mark5 Plus is manufactured under the controls of FDA-GMP's and ISO-9001 and EN46001 Quality Management Systems.

Anand Chitlangia  
President & CEO  
Nidek Medical Products, Inc.

This edition of the Mark5 Plus Maintenance Manual is effective for configuration models 1373-31100 and 1373-71100 beginning with serial number 992-18378. See Section 2 for configurations 1373-31103 and 1373-71103.

## Specifications (Stated at 70°F (20° C) and 1 Atm.)

|                                   |                       |
|-----------------------------------|-----------------------|
| Height.....                       | 26 in / 660 mm        |
| Width.....                        | 15 in / 380 mm        |
| Depth.....                        | 15 in / 380 mm        |
| Weight.....                       | 60 lbs / 27 kg        |
| Output.....                       | 90% +/- 3% @ 5 lpm    |
| Voltage.....                      | 115 V; 60 Hz          |
| Average Power Consumption.....    | 420W @ 60 Hz          |
| Cooling Air Temperature Rise..... | Less than 30°F / 16°C |

# Safety

## Oxygen Safety

Oxygen does not burn, but it does cause fire to burn more rapidly. Any material that will burn in air will ignite more rapidly and burn more violently in an oxygen-enriched atmosphere. While the oxygen concentrator does not contain a large volume of oxygen at any one given time, oxygen hazards can occur. If the output is allowed to return to the atmosphere no real hazard is likely; however, if the output hose is under a chair cushion or a pile of shop rags, these materials may become saturated with oxygen and could burn violently if ignited. The possibility of fire is the most serious hazard associated with the use of oxygen equipment. Observe the following precautions to prevent the possibility of a fire occurring in the vicinity of the oxygen outlet from this unit.

### CAUTION

#### **Eliminate Sources of Flame, Heat and Sparks**

A fire can be started by almost anything that produces flame, heat and sparks. Do not permit the use of matches, cigarette lighters, burning tobacco, candles, or other flame sources in rooms where oxygen is used or stored. Electrical heaters or equipment that gets hot or sparks during operation can be a source of ignition.

#### **Keep Flammable Materials Away From Oxygen**

It is not possible to completely eliminate the presence of combustible materials such as clothing, bedding, and even parts of the therapy equipment itself, but the presence of highly flammable materials can and must be avoided. Do not lubricate oxygen equipment. Oil and grease, including petroleum jelly such as Vaseline™ ignite very easily and burn vigorously in the presence of oxygen.

#### **Prevent Accumulation of Oxygen**

Since any material that will burn in air will burn more vigorously in an oxygen-enriched atmosphere, prevent any local accumulation of oxygen that may saturate combustibles. All connections on the apparatus should be made carefully and observed for leakage.

## Principles of Operation

The *Mark 5 Plus* Concentrator takes in room air that contains 21% oxygen, 78% nitrogen and trace amounts of inert gases. The air is filtered, compressed and passed through one of the two molecular sieve beds which allows the oxygen to pass but traps most of the nitrogen. Before the nitrogen saturates the molecular sieve in the first bed it is vented to the atmosphere and purged with oxygen while air is fed to the second bed. This process of alternate feeding and venting of beds is continued thus producing oxygen enriched gas at one end and exhaust at the other end of the sieve beds. The oxygen produced is then filtered, regulated, and set at a given flow rate for delivery to the patient.

The operation of the *Mark 5 Plus* can be grouped into three distinct systems; the *compressed air system*, the *process system* and the *delivery system*.

### Compressed Air System


Air enters the system through filters which prevent dust and other particles from entering and potentially damaging the air compressor. The compressor is a reciprocating two cylinder positive displacement type that compresses the air to the necessary pressure. The compressed air is then routed through a heat exchanger to reduce its temperature and then to the air surge tank. The surge tank reduces air pressure swing and traps moisture. The wick assembly, located at the bottom of the surge tank, allows the moisture to be evaporated to atmosphere. The system also has a safety relief valve to prevent pressure buildup should a blockage occur and a test port for measuring system pressure.

### Process System

The process system consist of two sieve beds and a rotary poppet valve. The valve controls the process by alternately feeding air into and exhausting nitrogen from the sieve beds. The valve utilizes a gear motor to cause a cam action to alternately open and close feed and exhaust poppets. Thus, what normally is a complex control system in other concentrators is reduced to a simple and easily understandable mechanical control in the *Mark 5 Plus*.

### Delivery System

The cyclic process described above causes pressure variations in the oxygen as it exits the beds. In order to provide a stable flow, the oxygen, after passing through a bacterial filter to remove any contaminants, is passed through a regulator to establish a fixed delivery pressure. A flow control valve consisting of multiple orifices then establishes a fixed flow rate to the patient. The use of a vacuum relief breaker and check valve in this section prevents the possibility of contaminants entering the system from the outside. Two pressure switches monitor the oxygen pressure and cause an alarm to sound if the pressure is high or low. An Oxygen Monitoring System for continuously monitoring oxygen purity is available as an option. The operation of the Oxygen Monitoring System is described separately in Section ② .



## Component Function

A brief description of the principal function of each of the major components is provided below:

|                          |  |
|--------------------------|--|
| CABINET FILTER.....      | Filters dust from cooling air                |
| INLET FILTER.....        | Filters Particles To Protect Compressor      |
| COMPRESSOR.....          | Pumps Air Into The System                    |
| RELIEF VALVE.....        | Safety Device To Limit System Pressur        |
| HEAT EXCHANGER.....      | Cools Air After Compression                  |
| WICK ASSY.....           | Removes Moisture from the Compressed Air     |
| AIR SURGE TANK.....      | Reduces air pressure fluctuations            |
| AIR TEST PORT.....       | A Connection For System Testing              |
| ROTARY POPPET VALVE..... | Controls The Pressure Swing Process          |
| SIEVE BEDS.....          | Separates Gases As Air Is Moved In And Out   |
| MUFFLER.....             | Reduces Exhaust Noise                        |
| OXYGEN TEST PORT.....    | A Connection For System Testing              |
| PRESSURE SWITCHES.....   | Activates Alarm For High Or Low Pressure     |
| BACTERIAL FILTER.....    | Removes contaminant from oxygen              |
| VACUUM BREAKER.....      | Prevents System Vacuum                       |
| REGULATOR.....           | Sets A Fixed Output Pressure                 |
| FLOW CONTROL VALV.E..... | Sets The Flow Rate To A Specific Flow        |
| CHECK VALVE.....         | Prevents Any Reverse Flow Into The System    |
| OUTLET.....              | Provides A Connection For Hose Or Humidifier |

# Operating Instructions

The operating instructions for the *Mark 5 Plus* Concentrator are provided in a separate booklet written for the patient's use. Maintenance personnel must make themselves familiar with the contents of this book and use it for all questions that relate to the operation of the unit.

## Initial Setup

Preparing a new *Mark 5 Plus* Concentrator for patient use consists of little more than taking it out of its shipping carton. After completing the Carton Unpacking Instructions found in the shipping carton follow these important checks to insure that shipping damage has not occurred and unit is performing properly.

### Step One: Flow Control Limit Adjustment

The flow control setting on the *Mark 5 Plus* Concentrator may be locked to limit the maximum output available to the user. This limit is preset at the factory to 5 liters per minute but may be changed to any flow setting. From the front of the unit remove the flow control knob by pulling the knob straight off its shaft. The back of the knob assembly has a flow locking disc. To reset the flow, remove the screw from the locking disc, lift and rotate the disc in the desired flow setting, reinstall the screw and replace the flow control knob.

### Step Two: Alarm Test

Before connecting power to the unit, press the ON button to test the battery and power failure alarm. It should sound when the switch is on and the power disconnected. If there is no alarm, check the 9-volt battery inside the front cabinet on switch plate and replace if necessary.

With the unit is plugged-in, push the ON button. The audible alarm should sound and the Service Required light should be lit and the compressor should run. After approximately 15 seconds both the alarm and the Serviced Required light should turn off.

### Step Three: Purity Test

With unit running connect an oxygen Analyzer (Example: OxyChek®) to oxygen outlet and verify results of 90% +/- 3% @ 5 LPM.

Mark 5 Plus Oxygen Concentrator is ready to be placed in service.

## Precautions and Safe Practices

The following precautions and safe practices are specifically directed to personnel who are experienced in the service of medical oxygen equipment. The device covered by these instructions is an electro-mechanical device that deserves all the respect inherent in servicing such equipment. Operation of this equipment also requires an awareness of all the potential hazards associated with handling and using oxygen.

## Maintenance Safety

Maintenance personnel must be properly trained in the service of oxygen concentrators. Electrical shock hazards exist inside the cabinet of this equipment. Good shop practices must be observed at all times.

**Use only approved replacement parts.** Medical oxygen equipment should be repaired with parts supplied by the manufacturer. Substitution of unauthorized parts can effect the performance of the unit. All parts that will come in contact with oxygen must be thoroughly cleaned. Remember, human health is affected by the end product of this device and a foreign substance left by maintenance could be inhaled by a patient. Thread sealant and lubricants used on oxygen equipment must be of approved quality for such use.

Always follow printed instructions when maintaining or repairing this equipment. Repair procedures may involve special techniques that must be followed when performing such work.

### **WARNING**

The sieve beds of the concentrator covered by these instructions contain a highly compressed spring. **DO NOT ATTEMPT TO DISASSEMBLE** a sieve bed or the spring may eject parts with dangerous velocity sufficient to cause personal injury.



## Preventive Maintenance

To ensure proper performance of the Mark 5 Plus Oxygen Concentrator purity checks are recommended every three (3) months.

When operating under normal conditions the following maintenance schedule is recommended and should be followed for proper operation of the unit. In severe conditions such as high temperature, dusty and/or dirty atmosphere, it may be necessary to perform the maintenance functions at more frequent intervals. Where procedures are involved, they are shown elsewhere in this publication.

### Preventive Maintenance Schedule Mark 5 Plus Oxygen Concentrator

|                  | 12 Months   | 24 Months | 36 Months |
|------------------|---|-----------|-----------|
| Inlet Filter     | X   | X         | X         |
| Bacterial Filter | X   | X         | X         |
| Battery          | X   | X         | X         |
| Muffler          |   |           | X         |
| Wick Assy        |   |           | x         |
| Compressor       | Rebuild or perform maintenance when oxygen purity and/or operating pressure falls below minimum specifications. |           |           |

*Note: In the Reference Material Section ③ of this manual we have provided a Maintenance Record Check List sample for your convenience.*

# Trouble Shooting


## System Problem Diagnosis

With system problem isolation in mind, the *Mark 5 Plus* was designed with two test ports that allow the monitoring of the system pressures at two key points. The pressure swings / ranges at these test ports will provide important diagnostic clues when used with the SYSTEM DIAGNOSIS chart. Each port has a plug that is held into place by a quick disconnect fitting. The plug may be removed by depressing the ring around the plug. A suitable pressure gauge with  $\frac{1}{4}$  diameter tube lead can then be connected to the quick disconnect fitting to measure the pressure swings.

### Normal Pressure Swings Are:

\*20 – 30 psig (140-205 kPa) at Air Test Port; 5 – 6 psig (35-42 kPa) difference between high pressure and low pressure

\*10 – 15 psig (70-105 kPa) at Oxygen Test Port; 1 – 2 psig (7-15 kPa) difference between high pressure and low pressure



Note: To obtain accurate readings please observe the following check list.


1. Make sure that the filters are clean.
2. Leak test all fittings and repair all leaks (a leak will distort all measurements).
3. Make sure that the oxygen flow is set correctly at 5 lpm and that actual flow is 5 lpm with an accurate flowmeter..
4. Always be sure of the accuracy of your test instruments.

## System Diagnosis Chart

| Symptom                           | Possible Cause  | Corrective Action  |
|-----------------------------------|---|--|
| Power Failure Alarm               | 1) No Electrical Power  | A) Check Power Source  |
| Service Req'd Light               | 1) Low Oxygen Pressure<br>2) High Oxygen Pressure                                   | A) See Trouble Shooting Guide  |
| Low Product Flow                  | 1) System Leak  | A) Test Product Pressure<br>B) Adjust Regulator                                      |
|                                   | 1) Clogged Bacterial Filter   | C) Change Inlet Filter   |
| High Product Flow                 | 1) Product Regulator  | A) Test Product Pressure<br>B) Adjust Regulator                                      |
| Noisy Compressor                  | 1) Worn Piston Seals<br>2) Rods or Bearings   | A) Replace Piston Seals or Replace Compressor  |
| Compressor Turns Off and Comes On | 1) Compressor Thermal Overload Switch Activated Due to Overheating                  | A) Check Fan and Cooling Airway Path<br>B) Replace Fan<br>C) Replace Compressor      |
| Fluctuating Flow                  | 1) Dirty or Defective Regulator<br>2) Leak on Oxygen Side<br>3) Low Oxygen Pressure | A) Check Operation<br>B) Clean or Replace Regulator<br>C) See Trouble Shooting Guide |

The trouble shooting guide which follows has been prepared to provide a logical sequence of measurements and evaluations that will lead to a probable conclusion as to the cause of a problem. As with all diagnostic efforts, the more highly trained and experienced individual is more likely to determine the root cause of any problem.

# Troubleshooting Guide




## Repair Procedures

Compressor Replacement  
(Compressor Model 2650CE37)

(Part # 7355-3535 – New)  
(Part # 7355-3537 – Rebuilt)

Carefully follow the step-by-step procedure described below:

1. Disconnect concentrator from electrical source and remove front cabinet enclosure (2 Phillips head screws in the handle recesses; 4 screws around the base).
2. Remove front of the compressor enclosure (6 Phillips head screws).
3. Disconnect hose from the exhaust muffler located on the compressor front enclosure.
4. Disconnect the capacitor wires from the capacitor located on the compressor front enclosure.
5. Remove the inlet hose from the right side of the compressor
6. Loosen top hose clamp on left side of the compressor and remove the hose from compressor outlet fitting.
7. Remove the two Phillips head screws securing the compressor plate.
8. Clip the tie wrap above the compressor box enclosure that bundles the compressor electrical wires.
9. Disconnect the compressor electrical wires from the terminal block located on the top of the compressor box enclosure.
10. Slide compressor forward approximately 6 mm, lift up and pull out.
11. Remove all fittings including the relief valve.
12. Remove the compressor plate from the compressor by removing the four bolts holding the compressor and compressor plate together.
13. Reinstall replacement compressor in reverse order using thread sealant on metal fittings.
14. Check for leaks.




## Sieve Module Assembly Replacement

(P/N 9012-8G06)

The Sieve Module Assembly in the *Mark 5 Plus* Concentrator consist of: 2 Sieve Beds, Gearmotor and Control Valve and is replaced by removing the entire assembly. Carefully follow the instructions below.

1. Disconnect the gearmotor wires from the terminal strip.
2. Disconnect product line tubing.
3. Disconnect exhaust hose from fitting.
4. Pull open Velcro™ straps from sieve beds.
5. Remove the air hose from the back of the control valve by loosening the clamp.
6. Install in reverse
7. Check for leaks.




### WARNING

The sieve beds of the concentrator covered by these instructions contain a highly compressed spring. DO NOT ATTEMPT TO DISASSEMBLE a sieve bed or the spring may eject parts with dangerous velocity sufficient to cause personal injury.


# Electrical Schematic

## Mark 5 Plus Oxygen Concentrator



# Parts Call-Out Illustration (Side View)


Mark 5 Plus Oxygen Concentrator





| MODEL        | ILL# | ITEM PART NUMBER | DESCRIPTION                             |
|--------------|------|------------------|---|
| 1373-31110T  | 1    | 5161-1027        | CAPACITOR 15 MFD                        |
| 1373-71100   | *    | 5161-1029        | CAPACITOR RUBBER CAP                    |
| 1373-71100C  | 2    | 5190-2233        | TIE WRAP HOLDER                         |
| 1373-71110   | *    | 5190-2234        | TIE WRAP 12 INCH LONG                   |
| 1373-71110CA | 3    | 6491-1005        | CONNECTOR, 4-WAY, 3/16 BARB             |
| 1373-71100JM | 4    | 6814-9228        | ELBOW NYLON DOUBLE BARB ½ INCH          |
| 1373-71110T  | 5    | 6816-4162        | TEE MALE BRANCH 1/8 NPT x ¼ ODT         |
| 1373-71110TB | 6    | 6814-4183        | ELBOW BRASS ¼ NPT x 3/8 HOSE BARB       |
| 1373-81110   | 7    | 6956-9673        | CHECK VALVE / VACUUM BREAKER            |
|              | 8    | 6998-9005        | FLOW CONTROL VALVE (1/8 to 5 LPM)       |
|              |      | 6998-9001        | FLOW CONTROL VALVE (1/2 to 6 LPM)       |
|              | 9    | 7206-0027        | BATTERY 9 VOLT                          |
|              | 10   | 7355-3535        | COMPRESSOR 115 V (NEW)                  |
|              | *    | 7355-3537        | COMPRESSOR 115 V (REBUILT)              |
|              | *    | 7355-3556        | COMPRESSOR REBUILD KIT                  |
|              | *    | 6953-9073        | RELIEF VALVE 44 PSIG                    |
|              | 11   | 7500-0004        | GEARMOTOR/VALVE ASSEMBLY                |
|              | *    | 7500-0003        | GEARMOTOR ASSEMBLY                      |
|              | 12   | 7603-0005A       | BLOWER ASSEMBLY                         |
|              | 13   | 7631-1053        | BACTERIAL FILTER                        |
|              | 14   | 7631-1058        | EXHAUST MUFFLER                         |
|              | 15   | 7631-1059        | INLET FILTER                            |
|              | 16   | 7631-1063        | CAP INLET FILTER                        |
|              | 17   | 7764-5587        | BUSHING ¾ INCH                          |
|              | 18   | 7854-6110        | HOSE BLUE PVC ¾ IN OD                   |
|              | 19   | 7854-6051        | HOSE SILICONE 5/32 ID x 3 LONG          |
|              | 20   | 7854-6049        | HOSE SILICONE 5/32 ID x 8 LONG          |
|              | 21   | 8300-8068        | CASTER                                  |
|              | 22   | 8316-5018        | HOURLY METER ASSEMBLY                   |
|              | 23   | 8354-3071        | VIBRATION ISOLATOR BLUE                 |
|              | 24   | 8816-1056        | REGULATOR                               |
|              | 25   | 8986-0016        | SWITCH ON/OFF                           |
|              |      | 8986-0023        | SWITCH ON/OFF W/10A BREAKER(SCHURTER)   |
|              | 26   | 9010-0G14        | SIEVE BED STRAP VELCRO                  |
|              | 27   | 9010-8G03        | SHOCK MOUNT KIT (COMPLETE)              |
|              | 28   | 9012-8G06        | SIEVE MODULE                            |
|              | *    | 9012-5G00        | SIEVE BED (NEW)                         |
|              | *    | 9012-5G01        | SIEVE BED (REBUILT)                     |
|              | 29   | 9012-7G00        | CONTROL VALVE ASSEMBLY                  |
|              | 30   | 0450-0001        | SWITCH PLATE ASSEMBLY (STANDARD SWITCH) |
|              |      | 0450-0011        | SWITCH PLATE ASSEMBLY (SCHURTER)        |
|              | *    | 9114-9G14        | SWITCH PLATE                            |
|              | *    | 8230-9037        | ALARM LIGHT                             |
|              | 31   | 9010-4G05        | AIR SURGE TANK                          |
|              | 32   | 9114-6G37        | CABINET BASE ASSEMBLY                   |
|              | 33   | 9114-6G20        | DUAL COIL HEAT EXCHANGER                |
|              | 34   | 9114-9G12        | RESONATOR                               |
|              | 35   | 9077-7008        | TUBING BLUE POLYURETHANE ¼ OD           |
|              | 36   | 6491-1004        | ELBOW 3/16 HOSE X 1/8 NPT NYLON         |
|              |      |                  |   |


# Parts Call-Out Illustration (Top View)



| MODEL        | ILL# | ITEM PART NUMBER | DESCRIPTION                     |
|--------------|------|------------------|---------------------------------|
| 1373-31110T  | 1    | 7603-0005A       | BLOWER ASSEMBLY (115 V)         |
| 1373-71100   | 2    | 5120-3181        | WIRE SADDLE WIDE                |
| 1373-71100C  | 3    | 6120-9855        | SCREW PHILLIPS HEAD BLACK       |
| 1373-71110   | 4    | 5161-1016        | TERMINAL STRIP                  |
| 1373-71110CA | 5    | 9010-4G01        | ELBOW/TUBE ASSEMBLY             |
| 1373-71100JM | 6    | 9010-4G05        | AIR SURGE TANK                  |
| 1373-71110T  | 7    | 7355-4775        | CONNECTOR 1/8 NPT x 1/4 ODT     |
| 1373-71110TB | 8    | 6814-9223        | ELBOW 3/8 NPT x 1/2 HOSE BARB   |
| 1373-81110   | 9    | 7854-6026        | HOSE 1/2 INCH BLACK x 20 LG     |
|              | 10   | 7342-1977        | CLAMP HOSE                      |
|              | 11   | 9007-4G25        | WICK ASSEMBLY                   |
|              | 12   | 7631-1063        | CAP INLET FILTER                |
|              | 13   | 7631-1059        | INLET FILTER                    |
|              | 14   | 8830-6140        | O-RING FOR FILTER HOLDER        |
|              | 15   | 9077-7008        | TUBING 1/4 OD BLUE POLYURETHANE |
|              | 16   | 6814-9228        | ELBOW 1/2 DOUBLE BARB NYLON     |
|              | 17   | 7854-6004        | EXHAUST HOSE 3/4 BLUE           |
|              | 18   | 6816-4162        | TEE 1/4 MALE BRANCH             |
|              | 19   | 8986-0012        | PRESSURE SWITCH LOW SILVER      |
|              | 20   | 8986-0015        | PRESSURE SWITCH HIGH GOLD       |
|              | 21   | 7854-6051        | HOSE 5/32 ID X 3 LONG           |
|              | 22   | 7854-6049        | HOSE 5/32 ID X 8 LONG           |
|              | 23   | 7854-6047        | HOSE 5/32 ID X 1.5 LONG         |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              |      |                  |                                 |
|              | *    |                  |                                 |
|              |      |                  |                                 |

Parts Call-Out Illustration (Front Outside View)

# Mark 5 Plus Oxygen Concentrator




| MODEL       | ILL# | ITEM PART NUMBER | DESCRIPTION               |
|-------------|------|------------------|---------------------------|
| 1373-31110T | 1    | 6120-9855        | SCREW PHILLIPS HEAD BLACK |

|              |    |           |                                      |
|--------------|----|-----------|--------------------------------------|
| 1373-71100   | 2  | 6132-4902 | SCREW PHILLIPS HEAD SILVER           |
| 1373-71100C  | 3  | 8300-8068 | CASTER                               |
| 1373-71110   | 4  | 9007-6G32 | FITTING HUMIDIFIER DISS 1240         |
| 1373-71110CA | 5  | 9014-6G10 | CABINET FRONT                        |
| 1373-71100JM | 6  | 9014-6G11 | CABINET BACK                         |
| 1373-71110T  | 7  | 9114-6G37 | CABINET BASE ASSEMBLY                |
| 1373-71110TB | 8  | 9114-6G45 | LABEL CABINET TOP                    |
| 1373-81110   | 9  | 9114-6G46 | LABEL ON / OFF                       |
|              | 10 | 9014-6G50 | FLOW CONTROL KNOB ASSY (1/8 – 5 LPM) |
|              |    | 9010-6G50 | FLOW CONTROL KNOB ASSY (1/2 – 6 LPM) |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |
|              |    |           |                                      |

Parts Call-Out Illustration (Back Outside View)

# Mark 5 Plus Oxygen Concentrator



| MODEL       | ILL# | ITEM PART NUMBER | DESCRIPTION               |
|-------------|------|------------------|---------------------------|
| 1373-31110T | 1    | 6120-9855        | SCREW PHILLIPS HEAD BLACK |

|              |    |            |                            |
|--------------|----|------------|----------------------------|
| 1373-71100   | 2  | 6132-4902  | SCREW PHILLIPS HEAD SILVER |
| 1373-71100C  | 3  | 7355-5019  | POWER CORD ASSEMBLY        |
| 1373-71110   | *  | 7281-5010  | STRAIN RELIEF              |
| 1373-71110CA | 4  | 8300-8068  | CASTER                     |
| 1373-71100JM | 5  | 9007-6G37  | CABINET FILTER             |
| 1373-71110T  | 6  | 9010-0G13  | CORD RETAINER              |
| 1373-71110TB | 7  | 9014-6G08  | ACCESS DOOR                |
| 1373-81110   | *  | 0580-0001  | LATCH ACCESS DOOR          |
|              | *  | 9114-6G52  | LABEL FILTER ACCESS DOOR   |
|              | 8  | 9014-6G10  | CABINET FRONT              |
|              | 9  | 9014-6G11  | CABINET BACK               |
|              | 10 | 9114-6G37  | CABINET BASE ASSEMBLY      |
|              | 11 | 9114-6G47  | BACK LABEL 115V 60 HZ      |
|              | *  | 9114-6G47C | BACK LABEL (CSA)           |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |
|              |    |            |                            |





# Operational Description for Oxygen Monitoring System

This portion of the manual is effective for configuration models 1373-31103 and 1373-71103 which are equipped with an Ultrasonic Oxygen Monitoring System. This is effective with for Oxygen Monitor equipped units beginning with serial number 98219700.

The Oxygen Monitor is an electronic module which continuously monitors the actual concentration of O<sub>2</sub> supplied by the **MARK 5 PLUS** concentrator. The Oxygen Monitor detects concentrations of 85% or below and triggers a visible and audible alarm. The Oxygen Monitor is available as an option. The Oxygen Monitor module reinforces and relays the standard operating alarms on the concentrator.

When powering on the **MARK 5 PLUS** the Oxygen Monitor module operates as below.

- In addition to the normal test of the **MARK 5 PLUS**, the Oxygen Monitor red indicator light comes on.
- The red indicator light may remain on for a few minutes (3 minutes maximum), which is the time necessary for the gas concentration supply to reach and exceed 85% O<sub>2</sub>.
- The green indicator light comes on after this value is reached and the concentrator is then operating satisfactorily.

## Operating principle

The oxygen monitor module is an electronic module that continuously monitors the oxygen enriched product that is actually supplied to the user.

The sensor, which utilizes the principle of sound velocity as a function of molecular weight, measures the velocity of sound in the product gas and compares it to stored values which correspond to the oxygen concentration. The preset threshold value is 85% for 115 volt devices.

A red light indicates a concentration below the preset threshold level. When the light is red for 15 minutes ( $\pm 2$  minutes) a continuously audible alarm is activated.

The oxygen monitor module incorporates a 9V transformer, a pressure sensor, the oxygen sensor, a 9V battery, red and green LEDs, the alarm buzzer and associated electronic circuits.

## Other alarms

Detection of absence of line voltage. In the event of a powerline cut, a continuous audible alarm is triggered (and the green LED goes out). The power supply for the alarm is a 9V battery.

Operating fault: In the case of a distribution malfunction (pressure of oxygen enriched air < 4 psig or > 20 psig) an audible alarm and visual alarm is triggered (red LED and continuous audible alarm).

## Interfacing

The oxygen monitoring system is fully contained on the printed circuit board on which the ON/OFF switch is located. Electrical power is supplied to the board from the switch terminals to terminals E1 and E2 on the board. The transformer and rectifier located on the board converts the line supplied power to the 5 V dc and 12 V dc power that is required to operate the electronic parts. Battery supplied power is received from the battery mounted on the board. This is the extent of the power requirement.

The product gas interface to the board is achieved by routing the product gas, which is to be delivered to the user via the outlet fitting, from the flow control valve to the oxygen monitoring unit. All of the gas produced by the concentrator is routed through the analyzer. From the analyzer sensor tube, the gas is routed to the concentrator outlet fitting and then to the user. The analyzer is designed to operate at low pressures only. The gas supplied must always pass through the regulator and the flow control valve before being routed to the analyzer tube. Use caution to be sure that no gas with pressure greater than 5 psig is fed to the analyzer sensor tube assembly.

Additional to the regulated product gas being connected to the sensor tube, a signal of unregulated product gas is connected to the pressure sensor located on the printed circuit board. This provides the signal to sound the low or high pressure alarm.

## Operation

The oxygen monitor is a single level oxygen monitoring system. The monitor is calibrated at both air and concentrator levels of oxygen concentration. The internal alarm switch is pre-set to a setting of 85 percent oxygen purity. When the concentration of the input gas drops below pre-set level, the red LED is lighted. If the concentration stays below the pre-set level for approximately 15 minutes, an audible alarm signal is generated. When the average purity is maintained above 85 percent for several concentrator cycles, the alarms are reset.


The process high and low pressure alarms are also included in the printed circuit board. It is no longer necessary to supply the high and low pressure switches in the product circuit.


## Trouble Shooting

- 1 Connect mains power to terminals E1 and E2.
- 2 Connect concentrator product gas to the left barb on the sensor tube. Check the purity of the product gas at the outlet fitting.
- 3 Turn power on. The alarm should sound until product pressure exceeds 4 psig.
- 4 Allow unit to reach operating temperature before starting the test. (1-3 minutes).
- 5 When gas purity reaches 85 percent, the green LED should light.
- 6 Bleed off some product gas at the sieve bed to reduce purity. When the input gas concentration is below the pre-set level, the red LED should be lighted. After about 15 minutes, the audible alarm should sound.

# Electrical Schematic

## Mark 5 Plus Oxygen Concentrator with Oxygen Monitoring System





Pneumatic outline diagram with Oxygen Monitor System

## Component Function

A brief description of the principal function of each of the major components is provided below:

|                          |  |
|--------------------------|--|
| CABINET FILTER.....      | Filters dust from cooling air              |
| INLET FILTER.....        | Filters Particles To Protect Compressor    |
| COMPRESSOR.....          | Pumps Air Into The System                  |
| RELIEF VALVE.....        | Safety Device To Limit System Pressure     |
| HEAT EXCHANGER.....      | Cools Air After Compression                |
| WICK ASSY.....           | Removes Moisture from the Compressed Air   |
| AIR SURGE TANK.....      | Reduces air pressure fluctuations          |
| AIR TEST PORT.....       | A Connection For System Testing            |
| ROTARY POPPET VALVE..... | Controls The Pressure Swing Process        |
| SIEVE BEDS.....          | Separates Gases As Air Is Moved In And Out |
| MUFFLER.....             | Reduces Exhaust Noise                      |
| OXYGEN TEST PORT.....    | A Connection For System Testing            |
| PRESSURE SENSOR.....     | Activates Alarm For High Or Low Pressure   |
| BACTERIAL FILTER.....    | Removes contaminant from oxygen            |
| VACUUM BREAKER.....      | Prevents System Vacuum                     |
| REGULATOR.....           | Sets A Fixed Output Pressure               |
| OXYGEN SENSOR.....       | Activates alarm for low purity             |
| FLOW CONTROL VALVE.....  | Sets The Flow Rate To A Specific Flow      |
| CHECK VALVE.....         | Prevents Any Reverse Flow Into The System  |
| OUTLET.....              | Provides Connection For Hose Or Humidifier |

## Preventive Maintenance

To ensure proper performance of the Mark 5 Plus Oxygen Concentrator purity checks are recommended every twelve (12) months or each time a unit is returned from a patient.

When operating under normal conditions the following maintenance schedule is recommended and should be followed for proper operation of the unit. In severe conditions such as high temperature, dusty and/or dirty atmosphere, it may be necessary to perform the maintenance functions at more frequent intervals. Where procedures are involved, they are shown elsewhere in this publication.

### Preventive Maintenance Schedule Mark 5 Plus Oxygen Concentrator with Oxygen Monitoring System


|                  | 12 Months   | 24 Months | 36 Months |
|------------------|---|-----------|-----------|
| Inlet Filter     | X   | x         | X         |
| Bacterial Filter | X   | x         | X         |
| Battery          | X   | x         | X         |
| Muffler          |   |           | X         |
| Wick             |   |           | X         |
| Compressor       | Rebuild or perform maintenance when oxygen purity and/or operating pressure falls below minimum specifications. |           |           |

The most frequent cause of low purity is that the air supply has been restricted. It is particularly important to be sure that the cabinet air filter is cleaned on a regular basis. If it is permitted to get dirty, it will reduce the amount of cooling air available to the machine which may result in premature loss of performance of the compressor and the sieve module.

*Note: In the Reference Material Section ③ of this manual we have provided a Maintenance Record Check List sample for your convenience.*

# Parts Call-Out Illustration (Side View)


Mark 5 Plus Oxygen Concentrator with Oxygen Monitoring System



| MODEL        | ILL# | ITEM PART NUMBER | DESCRIPTION                          |
|--------------|------|------------------|--------------------------------------|
| 1373-31110T  | 1    | 5161-1027        | CAPACITOR 15 MFD                     |
| 1373-71100   | *    | 5161-1029        | CAPACITOR RUBBER CAP                 |
| 1373-71100C  | 2    | 5190-2233        | TIE WRAP HOLDER                      |
| 1373-71110   | *    | 5190-2234        | TIE WRAP 12 INCH LONG                |
| 1373-71110CA | 3    | 6491-1005        | ADAPTER BARBED NYLON 4 WAY           |
| 1373-71100JM | 4    | 6814-9228        | ELBOW NYLON DOUBLE BARB ½ INCH       |
| 1373-71110T  | 5    | 6816-4162        | TEE MALE BRANCH 1/8 NPT x ¼ ODT      |
| 1373-71110TB | 6    | 6814-4183        | ELBOW BRASS ¼ NPT x 3/8 HOSE BARB    |
| 1373-81110   | 7    | 6956-9673        | CHECK VALVE / VACUUM BREAKER         |
|              | 8    | 6998-9005        | FLOW CONTROL VALVE                   |
|              | 9    | 7206-0027        | BATTERY 9 VOLT                       |
|              | 10   | 7355-3535        | COMPRESSOR 115 V (NEW)               |
|              | *    | 7335-3537        | COMPRESSOR 115 V (REBUILT)           |
|              | *    | 7355-3556        | COMPRESSOR REBUILD KIT )             |
|              | *    | 6953-9073        | RELIEF VALVE 44 PSIG                 |
|              | 11   | 7500-0004        | GEARMOTOR/VALVE ASSEMBLY             |
|              | *    | 7500-0003        | GEARMOTOR ASSEMBLY                   |
|              | 12   | 7603-0005A       | BLOWER ASSEMBLY                      |
|              | 13   | 7631-1053        | BACTERIAL FILTER                     |
|              | 14   | 7631-1058        | EXHAUST MUFFLER                      |
|              | 15   | 7631-1059        | INLET FILTER                         |
|              | 16   | 7631-1063        | CAP INLET FILTER                     |
|              | 17   | 7764-5587        | BUSHING ¾ INCH                       |
|              | 18   | 7854-6110        | HOSE BLUE PVC ¾ IN OD                |
|              | 19   | 7854-6049        | HOSE MILKY SILICONE 5/32 ID x 8 LONG |
|              | 20   | 7854-6047        | HOSE SILICONE 5/32 ID x 1.5 LONG     |
|              | 21   | 8300-8068        | CASTER                               |
|              | 22   | 8316-5018        | HOURMETER ASSEMBLY                   |
|              | 23   | 8354-3071        | VIBRATION ISOLATOR BLUE              |
|              | 24   | 8816-1066        | REGULATOR                            |
|              | 25   | 8986-0022        | SWITCH ON/OFF W/5A BREAKER(SCHURTER) |
|              | 26   | 9010-0G14        | SIEVE BED STRAP VELCRO               |
|              | 27   | 9010-8G03        | SHOCK MOUNT KIT (COMPLETE)           |
|              | 28   | 9012-8G06        | SIEVE MODULE ASSEMBLY                |
|              | *    | 9012-5G00        | SIEVE BED (NEW)                      |
|              | *    | 9012-5G01        | SIEVE BED (REBUILT)                  |
|              | 29   | 9012-7G00        | CONTROL VALVE ASSEMBLY               |
|              | 30   | 6500-4400        | OXYGEN MONITOR BOARD                 |
|              | 31   | 9010-4G05        | AIR SURGE TANK                       |
|              | 32   | 9114-6G37        | CABINET BASE ASSEMBLY                |
|              | 33   | 9114-6G20        | DUAL COIL HEAT EXCHANGER             |
|              | 34   | 9114-9G12        | RESONATOR                            |
|              | 35   | 9077-7008        | TUBING BLUE POLYURETHANE ¼ OD        |
|              |      |                  |                                      |

# Parts Call-Out Illustration (Top View)

Mark 5 Plus Oxygen Concentrator with Oxygen Monitoring System




| MODEL | ILL# | ITEM PART NUMBER | DESCRIPTION |
|-------|------|------------------|-------------|
|-------|------|------------------|-------------|



|              |    |            |   |
|--------------|----|------------|---|
| 1373-31110T  | 1  | 7603-0005A | BLOWER ASSEMBLY                         |
| 1373-71100   | 2  | 5120-3181  | WIRE SADDLE WIDE                        |
| 1373-71100C  | 3  | 6120-9855  | SCREW PHILLIPS HEAD BLACK               |
| 1373-71110   | 4  | 5161-1016  | TERMINAL STRIP                          |
| 1373-71110CA | 5  | 9010-4G01  | ELBOW/TUBE ASSEMBLY                     |
| 1373-71100JM | 6  | 9010-4G05  | AIR SURGE TANK                          |
| 1373-71110T  | 7  | 7355-4775  | CONNECTOR 1/8 NPT x 1/4 ODT             |
| 1373-71110TB | 8  | 6814-9223  | ELBOW 3/8 NPT x 1/2 HOSE BARB           |
| 1373-81110   | 9  | 7854-6026  | HOSE 1/2 INCH BLACK x 20 LG             |
|              | 10 | 7342-1977  | CLAMP HOSE                              |
|              | 11 | 9007-4G25  | WICK ASSEMBLY                           |
|              | 12 | 7631-1063  | CAP INLET FILTER                        |
|              | 13 | 7631-1059  | INLET FILTER                            |
|              | 14 | 8830-6140  | O-RING FOR FILTER HOLDER                |
|              | 15 | 9077-7008  | TUBING 1/4 OD BLUE POLYURETHANE         |
|              | 16 | 6814-9228  | ELBOW 1/2 DOUBLE BARB NYLON             |
|              | 17 | 7854-6004  | EXHAUST HOSE 3/4 BLUE                   |
|              | 18 | 6816-4162  | TEE 1/4 MALE BRANCH                     |
|              | 19 | 7854-6051  | HOSE 5/32 ID x 11/32 OD x 3 LG SILICONE |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              |    |            |   |
|              | *  |            |   |
|              |    |            |   |

Parts Call-Out Illustration (Front Outside View)  
 Mark 5 Plus Oxygen Concentrator with Oxygen Monitoring System



| MODEL | ILL# | ITEM PART NUMBER | DESCRIPTION |
|-------|------|------------------|-------------|
|-------|------|------------------|-------------|

